
HIPAP





User description on new functions and improvements for HiPAP/ μ PAP Transceiver software.

Last software version is HiPAP 4.0.2 -27. Mars 2019 for Win7/10

Release Description

HiPAP Release 4.0.2 - 27. Mars. 2019

New functions

• Enable reading of PSU voltage and temperature (x82)

Improvements

• Tested for Win10

Release History

HiPAP Release 3.5.4 - 17. Nov. 2017

New functions

NA

Improvements

• When no reply was present on one transponder in the LBL array, no position was computed and no position telegram sent to DP.

HiPAP Release 3.5.3 - 29. June 2017

New functions

• Added support for HiPAP 200

Improvements

NA.

HiPAP Release 3.5.2 - 27. June 2017

New functions

NA

Improvements

- Fixed a problem with receive filters for Remus channels on HiPAP 502.
- Improved startup when transceiver network is disconnected.

HiPAP Release 3.5.1 - 25. April 2017

New functions

NA

Improvements

• Fixed memory leak that could cause system to stop.

HiPAP Release 3.5.0 - 27. March 2017

New functions

- Improved range accuracy for Cymbal (M channels) and FSK (B channels.)
- Added support for FSK channel B33 and B44

Improvements

- Fixed a problem with receive filters for Remus channels on HiPAP 502.
- Improved robustness for cymbal transponder channels